



The all new Flexion from CJ-Optik Redefining excellence and economy, completely!!

Exclusive unibody construction offers purity of form and perfection of function. The microscope is created using a minimalist design philosophy, demanding that less is truly more. All essentials, Included!
The highest quality German optical and machining skills are combined to create a no-compromise system of elegance and simplicity.

Suddenly, all of the confusing options
become crystal clear!



Apo and Plan – apochromatic optics with the latest index-matched coatings are integrated, offering the highest correction for aberration, color and light transmission available. Precise and innovative construction provides superior ergonomic placement of all functions and controls, along with unparalleled ease-of-movement!



MonoGlobe, the weightless movement balancing system

Superior light transmission/throughput for optimal viewing and documentation



HD-Imaging ports are optimized for the newest digital cameras (APS-C/FullFrame)
Continuous power to all Sony Alpha cameras


CJ-Optik all new Flexion microscope!

Made in Germany,
with precision.

All of the essentials,
optimized in every way
for the best possible
optical, ergonomic and
mechanical performance,
just as you demand!



Elegance, high-performance,
solid build and simplicity, all at once



USB, HDMI, Medical monitor power
integrated for a better workflow

Ergotube 30° with binocular inclination
0-200° and long eye relief eyepieces